

Manual



Tierco[®]
by Seating Matters

Manual Tilt in Space

Toileting and Showering Chair



Free Assessment

www.tierco.com
info@tierco.com

Manual Tilt in Space

Toileting and Shower Chair

The Manual Tilt in Space toileting and showering chair is designed to fit over Clos-o-mat and Geberit toilets plus 90% of all domestic toilets. The Tilt in Space is a product that truly has both the client and carers needs firmly in mind. To appreciate the flexibility of this product we would welcome the opportunity to demonstrate this product to you in person. Designed by Daily Care Ltd, all products are manufactured in the UK to the highest standards.

Adjustable: Reclines using a break leaver to allow the chair to tilt to 30°


Manoeuvrable: Easily transported from room to room.

Detachable: Foot plates that can be individually adjusted and detached. Side arms are also detachable for side transfers.

Mesh: Easily removable so that it can be washed. The mesh is coated with Polygiene to avoid cross contamination.

Optional Accessories: Pads, rolls, bib harness, aperture, head supports, headrest and commode facility.

Construction: All stainless steel, can be used in showers, bathrooms and wet rooms.



Dimensions	
A: Headrest Width	298mm
B: Back Height	635mm
C: Thigh Length	492mm
D: Seat Height	530mm
E: Seat Width	525mm
F: Base Width	645.5mm

The various dimensions shown in the diagram are the number of flexible measurements which can be modified to suit the user needs.

30° Tilt Angle



Tierco[®]
by Seating Matters